


## (2S,4R)-Ethyl 5-([1,1'-biphenyl]-4-yl)-4-amino-2-methylpentanoate hydrochloride

## Chemical Properties

CAS No. :	2259707-96-3
Formula:	C <sub>20</sub> H <sub>26</sub> ClNO <sub>2</sub>
Molecular Weight:	347.88
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(2S,4R)-Ethyl 5-([1,1'-biphenyl]-4-yl)-4-amino-2-methylpentanoate hydrochloride is an intermediate in the synthesis of various active compounds.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8746 mL	14.3728 mL	28.7455 mL
5 mM	0.5749 mL	2.8746 mL	5.7491 mL
10 mM	0.2875 mL	1.4373 mL	2.8746 mL
50 mM	0.0575 mL	0.2875 mL	0.5749 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481