


## Oxiran-2-ylmethyl 4-nitrobenzenesulfonate

## Chemical Properties

CAS No. :	158505-66-9
Formula:	C <sub>9</sub> H <sub>9</sub> NO <sub>6</sub> S
Molecular Weight:	259.24
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Oxiran-2-ylmethyl 4-nitrobenzenesulfonate is an intermediate used in the synthesis of various active compounds.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8574 mL	19.2871 mL	38.5743 mL
5 mM	0.7715 mL	3.8574 mL	7.7149 mL
10 mM	0.3857 mL	1.9287 mL	3.8574 mL
50 mM	0.0771 mL	0.3857 mL	0.7715 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481