

JS006

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	JS006 is a humanized IgG4k monoclonal antibody inhibitor targeting TIGIT/CD155. It is applicable in the study of advanced and/or metastatic cancers, such as non-small cell lung cancer (NSCLC), advanced lymphoma, and advanced/metastatic triple-negative breast cancer (TNBC), particularly when resistance to PD-1/PD-L1 inhibitors has developed.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481