

S-095029

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	S-095029 is a humanized IgG1 monoclonal antibody inhibitor that targets NKG2A. It effectively reduces Fc effector functions and inhibits interaction with its ligand, HLA-E, while enhancing the antibody-dependent cellular cytotoxicity mediated by other intact Fc function antibodies. S-095029 demonstrates potent antitumor activity, increasing the cytotoxic and cytokine (IFN $\gamma$ , TNF- $\alpha$ , and CXCL9) secretion capabilities of NK cells and $\gamma\delta$ T cells when co-cultured with cancer cells.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481