

HLX26

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	HLX26 is a humanized anti-LAG-3 monoclonal antibody with antitumor properties, making it particularly useful for advanced solid tumors. Additionally, HLX26 can be applied in lymphoma research.
Targets(IC50)	LAG-3

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481