


(Rac)-IBR120

Chemical Properties

CAS No. : 1996673-44-9
 Formula: C₂₃H₂₀N₂O₂S
 Molecular Weight: 388.48
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	(Rac)-IBR120 is the racemate of IBR120, a RAD51 inhibitor that significantly suppresses the growth of MDA-MB-468 cells. It is applicable in research related to treatment-resistant cancers, such as triple-negative breast cancer.
Targets(IC50)	DNA/RNA Synthesis

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5741 mL	12.8707 mL	25.7414 mL
5 mM	0.5148 mL	2.5741 mL	5.1483 mL
10 mM	0.2574 mL	1.2871 mL	2.5741 mL
50 mM	0.0515 mL	0.2574 mL	0.5148 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481