


3-(3,5-Dichlorophenyl)-1-methylhydantoin

Chemical Properties

CAS No. :	27387-90-2
Formula:	C ₁₀ H ₈ Cl ₂ N ₂ O ₂
Molecular Weight:	259.09
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-(3,5-Dichlorophenyl)-1-methylhydantoin is a pharmaceutical intermediate utilized in the synthesis of various active compounds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8597 mL	19.2983 mL	38.5966 mL
5 mM	0.7719 mL	3.8597 mL	7.7193 mL
10 mM	0.386 mL	1.9298 mL	3.8597 mL
50 mM	0.0772 mL	0.386 mL	0.7719 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481