

Anti-CD20 Antibody (mAb 1.5.3)

Chemical Properties

CAS No. : 922534-56-3

Formula:

Molecular Weight:

Store at low temperature

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-CD20 Antibody (mAb 1.5.3) is a fully human IgG1 anti-CD20 antibody. It enhances pro-apoptotic activity in vitro and mediates complement-dependent and antibody-dependent cytotoxic effects. In various tumor xenograft models, this antibody shows increased anti-tumor activity. In primate pharmacodynamics models, it exhibits excellent B cell depletion in lymphatic organs and bone marrow. It is applicable in the study of B cell malignancies.
Targets(IC50)	Apoptosis, Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481