

SFN68

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SFN68 is a humanized IgG1 neutralizing monoclonal antibody targeting HGF. It significantly inhibits the interaction between HGF and c-Met, the activation of c-Met, and HGF-mediated cell scattering and proliferation. SFN68 is applicable for cancer research.
Targets(IC50)	c-Met/HGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481