


(2S,5S)-Methyl 5-phenylpyrrolidine-2-carboxylate hydrochloride

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 2170170-16-6 |
| Formula: | C ₁₂ H ₁₆ ClNO ₂ |
| Molecular Weight: | 241.71 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA. |



Biological Description

| | |
|---------------|--|
| Description | (2S,5S)-Methyl 5-phenylpyrrolidine-2-carboxylate hydrochloride serves as a pharmaceutical intermediate employed in synthesizing a variety of active compounds. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 4.1372 mL | 20.6859 mL | 41.3719 mL |
| 5 mM | 0.8274 mL | 4.1372 mL | 8.2744 mL |
| 10 mM | 0.4137 mL | 2.0686 mL | 4.1372 mL |
| 50 mM | 0.0827 mL | 0.4137 mL | 0.8274 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481