

MK-2060

Chemical Properties

CAS No. : 2812480-88-7

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	MK-2060 is a humanized monoclonal antibody inhibitor targeting FXI. It is applicable for research in atherosclerotic thrombosis and venous thromboembolism.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481