


Ethyl 4-(2-chloropyridin-3-ylamino)piperidine-1-carboxylate

Chemical Properties

CAS No. :	906371-78-6
Formula:	C ₁₃ H ₁₈ ClN ₃ O ₂
Molecular Weight:	283.75
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ethyl 4-(2-chloropyridin-3-ylamino)piperidine-1-carboxylate [1-Piperidinecarboxylic acid, 4-[(2-chloro-3-pyridinyl)amino]-, ethyl ester] is a pharmaceutical intermediate utilized in the synthesis of various active compounds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5242 mL	17.6211 mL	35.2423 mL
5 mM	0.7048 mL	3.5242 mL	7.0485 mL
10 mM	0.3524 mL	1.7621 mL	3.5242 mL
50 mM	0.0705 mL	0.3524 mL	0.7048 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481