


## 3-(4,4,5,5-Tetramethyl-1,3,2-dioxaborolan-2-yl)benzenesulfonyl fluoride

## Chemical Properties

CAS No. :	2393030-89-0
Formula:	C <sub>12</sub> H <sub>16</sub> BF <sub>4</sub> O <sub>4</sub> S
Molecular Weight:	286.13
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-(4,4,5,5-Tetramethyl-1,3,2-dioxaborolan-2-yl)benzenesulfonyl fluoride is a pharmaceutical intermediate utilized in the synthesis of a variety of active compounds.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4949 mL	17.4746 mL	34.9491 mL
5 mM	0.699 mL	3.4949 mL	6.9898 mL
10 mM	0.3495 mL	1.7475 mL	3.4949 mL
50 mM	0.0699 mL	0.3495 mL	0.699 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481