

BLX-301

## Chemical Properties

CAS No. : 2170268-76-3

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	BLX-301 is a humanized monoclonal antibody against CD20, exhibiting antitumor activity. It is applicable in studies related to diseases such as non-Hodgkin's lymphoma and rheumatoid arthritis.
Targets(IC50)	Immunology/Inflammation related

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481