


## 3-(4-Chlorophenoxy)-4-methylaniline

## Chemical Properties

CAS No. :	1400872-19-6
Formula:	C <sub>13</sub> H <sub>12</sub> ClNO
Molecular Weight:	233.69
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-(4-Chlorophenoxy)-4-methylaniline is a pharmaceutical intermediate used in the synthesis of various active compounds.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2792 mL	21.3959 mL	42.7917 mL
5 mM	0.8558 mL	4.2792 mL	8.5583 mL
10 mM	0.4279 mL	2.1396 mL	4.2792 mL
50 mM	0.0856 mL	0.4279 mL	0.8558 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481