

6-Chloropyrimidine-4-carboxylic acid

Chemical Properties

CAS No. :	37131-91-2
Formula:	C ₅ H ₃ ClN ₂ O ₂
Molecular Weight:	158.543
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	6-Chloropyrimidine-4-carboxylic acid is a pharmaceutical intermediate used in the synthesis of various active compounds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.3076 mL	31.5378 mL	63.0756 mL
5 mM	1.2615 mL	6.3076 mL	12.6151 mL
10 mM	0.6308 mL	3.1538 mL	6.3076 mL
50 mM	0.1262 mL	0.6308 mL	1.2615 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481