

mPEG14-Alcohol

Chemical Properties

CAS No. :	1059605-02-5
Formula:	C ₂₉ H ₆₀ O ₁₅
Molecular Weight:	648.78
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	mPEG14-Alcohol (mPEG14-OH) is a PEG derivative featuring a hydroxyl (OH) functional group, commonly utilized in drug delivery applications.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5414 mL	7.7068 mL	15.4135 mL
5 mM	0.3083 mL	1.5414 mL	3.0827 mL
10 mM	0.1541 mL	0.7707 mL	1.5414 mL
50 mM	0.0308 mL	0.1541 mL	0.3083 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481