

2-Amino-4-bromo-3-nitropyridine

Chemical Properties

CAS No. :	84487-10-5
Formula:	C ₅ H ₄ BrN ₃ O ₂
Molecular Weight:	218.01
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Amino-4-bromo-3-nitropyridine is an intermediate used in pharmaceutical synthesis to produce a variety of active compounds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5869 mL	22.9347 mL	45.8695 mL
5 mM	0.9174 mL	4.5869 mL	9.1739 mL
10 mM	0.4587 mL	2.2935 mL	4.5869 mL
50 mM	0.0917 mL	0.4587 mL	0.9174 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481