

Methyl (3S,6R)-6-methyl-3-piperidinecarboxylate L-tartaric acid salt

Chemical Properties


CAS No. : 1227911-35-4

Formula: C₁₂H₂₁N₀O₈

Molecular Weight: 307.30

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Methyl (3S,6R)-6-methyl-3-piperidinecarboxylate L-tartaric acid salt is a pharmaceutical intermediate used in the synthesis of various active compounds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2541 mL	16.2707 mL	32.5415 mL
5 mM	0.6508 mL	3.2541 mL	6.5083 mL
10 mM	0.3254 mL	1.6271 mL	3.2541 mL
50 mM	0.0651 mL	0.3254 mL	0.6508 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481