


2-Acetamidomalonnate sodium

Chemical Properties

CAS No. :	117976-12-2
Formula:	C ₅ H ₅ NNa ₂ O ₅
Molecular Weight:	205.08
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Acetamidomalonnate sodium is an intermediate used in the synthesis of various active compounds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8761 mL	24.3807 mL	48.7615 mL
5 mM	0.9752 mL	4.8761 mL	9.7523 mL
10 mM	0.4876 mL	2.4381 mL	4.8761 mL
50 mM	0.0975 mL	0.4876 mL	0.9752 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481