


Sonvuterkib

Chemical Properties

CAS No. :	2890225-50-8
Formula:	C ₂₃ H ₂₂ N ₈ O ₂ S
Molecular Weight:	474.54
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sonvuterkib (WX001) is a potent and orally active inhibitor of extracellular signal-regulated kinases (ERK), exhibiting IC ₅₀ values of 1.4 nM for ERK1 and 0.54 nM for ERK2. It effectively inhibits cell proliferation and possesses anticancer properties.
Targets(IC ₅₀)	ERK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1073 mL	10.5365 mL	21.073 mL
5 mM	0.4215 mL	2.1073 mL	4.2146 mL
10 mM	0.2107 mL	1.0537 mL	2.1073 mL
50 mM	0.0421 mL	0.2107 mL	0.4215 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481