

## Allyl-PEG3-OH

## Chemical Properties


CAS No. : 26150-05-0

Formula: C<sub>9</sub>H<sub>18</sub>O<sub>4</sub>

Molecular Weight: 190.24

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Allyl-PEG3-OH (Acryloyl-PEG3-OH) is a PEG derivative featuring a hydroxyl (OH) functional group, utilized in drug delivery applications.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2565 mL	26.2826 mL	52.5652 mL
5 mM	1.0513 mL	5.2565 mL	10.513 mL
10 mM	0.5257 mL	2.6283 mL	5.2565 mL
50 mM	0.1051 mL	0.5257 mL	1.0513 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481