


Diethyl 2-bromoglutarate

Chemical Properties

CAS No. : 7209-00-9
 Formula: C₉H₁₅BrO₄
 Molecular Weight: 267.12
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Diethyl 2-bromoglutarate is a pharmaceutical intermediate used in the synthesis of various active compounds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7436 mL	18.7182 mL	37.4364 mL
5 mM	0.7487 mL	3.7436 mL	7.4873 mL
10 mM	0.3744 mL	1.8718 mL	3.7436 mL
50 mM	0.0749 mL	0.3744 mL	0.7487 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481