


2-(DiMethylphosphoryl)benzenamine

Chemical Properties

CAS No. :	1197953-47-1
Formula:	C ₈ H ₁₂ NOP
Molecular Weight:	169.16
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-(DiMethylphosphoryl)benzenamine is a chemical intermediate used in the synthesis of various active compounds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.9116 mL	29.5578 mL	59.1156 mL
5 mM	1.1823 mL	5.9116 mL	11.8231 mL
10 mM	0.5912 mL	2.9558 mL	5.9116 mL
50 mM	0.1182 mL	0.5912 mL	1.1823 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481