


o-Tolidine Sulfone(contains 2,6-Dimethyl isomer)

Chemical Properties

CAS No. :	55011-44-4
Formula:	C ₁₄ H ₁₄ N ₂ O ₂ S
Molecular Weight:	274.34
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	o-Tolidine Sulfone (contains 2,6-Dimethyl isomer) [3,7-Diamino-2,8-dimethyldibenzothiophene sulfone (contains 2,6-Dimethyl isomer)] acts as a pharmaceutical intermediate in the synthesis of various active compounds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6451 mL	18.2256 mL	36.4511 mL
5 mM	0.729 mL	3.6451 mL	7.2902 mL
10 mM	0.3645 mL	1.8226 mL	3.6451 mL
50 mM	0.0729 mL	0.3645 mL	0.729 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481