

KHK6640

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	KHK6640 is a humanized antibody specific to amyloid beta oligomers. It demonstrates high efficacy in improving cognitive function in various rodent models of Alzheimer's disease. KHK6640 is applicable for research on neurological disorders such as Alzheimer's disease.
Targets(IC50)	Beta Amyloid

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481