


## CoA-S-acetyl-5-bromotryptamine

## Chemical Properties

CAS No. :	252860-65-4
Formula:	C <sub>33</sub> H <sub>47</sub> BrN <sub>9</sub> O <sub>17</sub> P <sub>3</sub> S
Molecular Weight:	1046.67
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	CoA-S-acetyl-5-bromotryptamine (Coenzyme A-S-acetyl 5-bromotryptamine) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9554 mL	4.7771 mL	9.5541 mL
5 mM	0.1911 mL	0.9554 mL	1.9108 mL
10 mM	0.0955 mL	0.4777 mL	0.9554 mL
50 mM	0.0191 mL	0.0955 mL	0.1911 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481