

TQB2916

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	TQB2916 is a humanized IgG2 monoclonal antibody agonist that targets CD40. It exhibits significant antitumor activity by occupying CD40 and activating immune function. TQB2916 is applicable for research in advanced solid tumors and lymphomas.
Targets(IC50)	TNF

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481