

AMG-167

## Chemical Properties

CAS No. : 1798349-30-0

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	AMG-167 is a monoclonal antibody targeting sclerostin, utilized in cancer research. Additionally, AMG-167 can be applied in studies related to osteoporosis.
Targets(IC50)	Wnt/beta-catenin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481