


3-Oxoicosanoyl-CoA

Chemical Properties

CAS No. :	88297-33-0
Formula:	C41H72N7O18P3S
Molecular Weight:	1076.03
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-Oxoicosanoyl-CoA (3-Oxoicosanoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9293 mL	4.6467 mL	9.2934 mL
5 mM	0.1859 mL	0.9293 mL	1.8587 mL
10 mM	0.0929 mL	0.4647 mL	0.9293 mL
50 mM	0.0186 mL	0.0929 mL	0.1859 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481