


Axareotide

Chemical Properties

CAS No. :	2126833-17-6
Formula:	C ₅₄ H ₆₈ ClN ₁₁ O ₁₂ S ₂
Molecular Weight:	1162.77
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Axareotide is an analog of somatostatin (Somatostatin).
Targets(IC50)	Somatostatin

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.860 mL	4.3001 mL	8.6002 mL
5 mM	0.172 mL	0.860 mL	1.720 mL
10 mM	0.086 mL	0.430 mL	0.860 mL
50 mM	0.0172 mL	0.086 mL	0.172 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481