


## Adipoyl-CoA

## Chemical Properties

CAS No. : 25119-42-0  
 Formula: C<sub>27</sub>H<sub>44</sub>N<sub>7</sub>O<sub>19</sub>P<sub>3</sub>S  
 Molecular Weight: 895.66  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Adipoyl-CoA (Adipoyl-Coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1165 mL	5.5825 mL	11.165 mL
5 mM	0.2233 mL	1.1165 mL	2.233 mL
10 mM	0.1116 mL	0.5582 mL	1.1165 mL
50 mM	0.0223 mL	0.1116 mL	0.2233 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481