

## cis-2-(Methoxycarbonyl)cyclopentanecarboxylic acid

## Chemical Properties

CAS No. :	92541-43-0
Formula:	C <sub>8</sub> H <sub>12</sub> O <sub>4</sub>
Molecular Weight:	172.179
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	cis-2-(Methoxycarbonyl)cyclopentanecarboxylic acid is an intermediate used in the synthesis of various active compounds.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8079 mL	29.0394 mL	58.0788 mL
5 mM	1.1616 mL	5.8079 mL	11.6158 mL
10 mM	0.5808 mL	2.9039 mL	5.8079 mL
50 mM	0.1162 mL	0.5808 mL	1.1616 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481