

MEDI-5884

Chemical Properties

CAS No. : 2643310-75-0

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	MEDI-5884 is a humanized IgG4Pk neutralizing monoclonal antibody that targets endothelial lipase (EL). By inhibiting EL, MEDI-5884 leads to an increase in high-density lipoprotein (HDL) levels, particularly phosphatidylinositol and cholesterol ester, thereby enhancing its functionality. This compound is applicable in research related to cardiovascular diseases, such as coronary artery disease (CAD).
Targets(IC50)	Lipase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481