

CTM01

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	CTM01 is a murine-derived IgG1 (or humanized IgG4) monoclonal antibody inhibitor targeting MUC1. It exhibits broad-spectrum anticancer activity against epithelial-derived solid tumors such as breast cancer, lung cancer, and ovarian cancer.
Targets(IC50)	Mucin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481