


Me-diazirine-cooh

Chemical Properties

CAS No. :	25055-86-1
Formula:	C ₅ H ₈ N ₂ O ₂
Molecular Weight:	128.13
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Me-diazirine-cooh [3-(3-Methyl-3H-diazirin-3-yl)propanoic acid] is a pharmaceutical intermediate used in the synthesis of various active compounds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.8046 mL	39.0229 mL	78.0457 mL
5 mM	1.5609 mL	7.8046 mL	15.6091 mL
10 mM	0.7805 mL	3.9023 mL	7.8046 mL
50 mM	0.1561 mL	0.7805 mL	1.5609 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481