


## 2-Chloro-5-(trifluoromethyl)pyrimidine

## Chemical Properties

CAS No. :	69034-12-4
Formula:	C <sub>5</sub> H <sub>2</sub> ClF <sub>3</sub> N <sub>2</sub>
Molecular Weight:	182.53
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-Chloro-5-(trifluoromethyl)pyrimidine is a pharmaceutical intermediate used in the synthesis of various active compounds.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.4786 mL	27.3928 mL	54.7855 mL
5 mM	1.0957 mL	5.4786 mL	10.9571 mL
10 mM	0.5479 mL	2.7393 mL	5.4786 mL
50 mM	0.1096 mL	0.5479 mL	1.0957 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481