

ABT-147

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ABT-147 is a humanized IgG1 monoclonal antibody inhibitor targeting IL-12. It is utilized in researching autoimmune diseases such as rheumatoid arthritis (RA), inflammatory diseases like inflammatory bowel disease (IBD), and other conditions including asthma.
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481