

4-Oxo-octanoyl-CoA

Chemical Properties


CAS No. : 2245809-15-6

Formula: C₂₉H₄₈N₇O₁₈P₃S

Molecular Weight: 907.71

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-Oxo-octanoyl-CoA (4-Oxo-octanoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1017 mL	5.5084 mL	11.0167 mL
5 mM	0.2203 mL	1.1017 mL	2.2033 mL
10 mM	0.1102 mL	0.5508 mL	1.1017 mL
50 mM	0.022 mL	0.1102 mL	0.2203 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481