

S095024

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	S095024 is a humanized IgG1 monoclonal antibody inhibitor targeting CD73. It is applicable in research related to acute respiratory distress syndrome (ARDS), non-small cell lung cancer (NSCLC), and advanced solid tumors [malignant neoplasms].
Targets(IC50)	CD73

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481