


## 1-Methyl-4-(2,4,6-trimethoxyphenyl)piperidin-3-ol

## Chemical Properties

|                   |   |
|-------------------|---|
| CAS No. :         | 234771-34-7   |
| Formula:          | C <sub>15</sub> H <sub>23</sub> NO <sub>4</sub>   |
| Molecular Weight: | 281.35  |
| Storage:          | Powder: -20°C for 3 years   In solvent: -80°C for 1 year<br>Actual storage temperature shall be subject to the COA. |



## Biological Description

|               |   |
|---------------|---|
| Description   | 1-Methyl-4-(2,4,6-trimethoxyphenyl)piperidin-3-ol is a pharmaceutical intermediate utilized in the synthesis of various active compounds. |
| Targets(IC50) | Others  |

## Preparing Stock Solutions

|       | 1mg       | 5mg        | 10mg       |
|-------|-----------|------------|------------|
| 1 mM  | 3.5543 mL | 17.7715 mL | 35.5429 mL |
| 5 mM  | 0.7109 mL | 3.5543 mL  | 7.1086 mL  |
| 10 mM | 0.3554 mL | 1.7771 mL  | 3.5543 mL  |
| 50 mM | 0.0711 mL | 0.3554 mL  | 0.7109 mL  |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481