


## 4-Amino-3-nitropyridin-2-ol

## Chemical Properties

CAS No. :	88511-57-3
Formula:	C <sub>5</sub> H <sub>5</sub> N <sub>3</sub> O <sub>3</sub>
Molecular Weight:	155.11
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4-Amino-3-nitropyridin-2-ol is a pharmaceutical intermediate employed in the synthesis of various active compounds.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.447 mL	32.2352 mL	64.4704 mL
5 mM	1.2894 mL	6.447 mL	12.8941 mL
10 mM	0.6447 mL	3.2235 mL	6.447 mL
50 mM	0.1289 mL	0.6447 mL	1.2894 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481