

Anti-EGFR Antibody (D2C7)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-EGFR Antibody (D2C7) is a monoclonal antibody targeting EGFR. The recombinant immunotoxin D2C7-(scdsFv)-PE38KDEL, developed from Anti-EGFR Antibody (D2C7), exhibits antitumor activity in glioblastoma models.
Targets(IC50)	EGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481