

Mocertatug

Chemical Properties

CAS No. : 3042847-92-4

Formula:

Molecular Weight:

Storage: Store at low temperature

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Mocertatug is an IgG1 κ -type humanized antibody that targets B7-H4/VTCN1, with its isotype control being HumanIgG1kappa, Isotype Control.
Targets(IC50)	Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481