

Velvakiment

Chemical Properties

CAS No. : 3036636-75-3

Formula:

Molecular Weight:

Storage: Store at low temperature

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Velvakiment is a humanized VHH nanobody that targets IL-17A. It is applicable in research on autoimmune diseases such as multiple sclerosis and rheumatoid arthritis.
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481