


BCL6 ligand-4

Chemical Properties

CAS No. :	2257479-50-6
Formula:	C ₂₁ H ₂₅ ClN ₆ O ₃
Molecular Weight:	444.92
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	BCL6ligand-4 is a BCL6 ligand used in the synthesis of PROTACs, such as BCL6PROTAC 1.
Targets(IC50)	Bcl-6,Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2476 mL	11.238 mL	22.476 mL
5 mM	0.4495 mL	2.2476 mL	4.4952 mL
10 mM	0.2248 mL	1.1238 mL	2.2476 mL
50 mM	0.045 mL	0.2248 mL	0.4495 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481