

Zarutatug

Chemical Properties

CAS No. : 3038456-02-6

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Zarutatug is a humanized IgG1 κ antibody with antitumor activity targeting CDH17 (cadherin 17). The recommended isotype control for this antibody is Human IgG1kappa, Isotype Control.
Targets(IC50)	Cadherin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481