

Vasigrobart

Chemical Properties

CAS No. : 3030730-59-4

Formula:

Molecular Weight:

Store at low temperature

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Vasigrobart is a humanized antibody expressed in CHO cells, targeting PDGFB. It consists of a hulgG1 heavy chain and a huk light chain, with a predicted molecular weight (MW) of 145 kDa. The isotype control for Vasigrobart is HumanIgG1kappa, Isotype Control. |
| Targets(IC50) | PDGFR |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481