

Gamgertamig

Chemical Properties

CAS No. : 3038119-48-8

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Gamgertamig is a humanized bispecific antibody of the IgG4 type that targets TNFRSF17 and CD3E.
Targets(IC50)	Integrin,TNF

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481