

MDX1201

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	MDX1201 is a humanized monoclonal antibody targeting PSMA1/7. Conjugated with the fluorescent dye AlexaTM488, MDX1201 serves as a probe for near-infrared fluorescence imaging (NIRF) of prostate cancer.
Targets(IC50)	PSMA

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481